

Auto Inflation Blood Pressure Monitor

UA-782



AND

A&D Medical

<http://www.aandd.jp>

Technology You Can Trust

Blood Pressure Monitor

Designed for Women

Model UA-782

Smooth and fast measurements

You may feel a little uncomfortable during blood pressure measurements because of the pressure applied to the arm and the noise while inflating the cuff. With our newest technology, those discomforts are a thing of the past. The UA-782 requires a shorter measurement time and its newly developed semi-hard shell cuff also caters to a wider range of arm sizes. The new whisper-silent pump further aids smooth and fast measurements. Arm range: 17cm-32cm

Irregular Heart Beat (IHB) indicator

This indicator illuminates if an IHB is detected during a measurement. The indicator alerts you of an irregular heartbeat, sometimes a symptom of cardiovascular disease. (Note: This is for warning purposes only. A&D recommends contacting a physician if the user observes this indicator frequently).

Blood pressure classification indicator

The World Health Organization (WHO) distinguishes 6 different blood pressure stages. This monitor automatically indicates your classification (from Optimal to Severe hypertension) after each measurement is completed.

90-reading memory

The UA-782 has a 90-reading memory capacity for monitoring your measurement trends. An average of all the stored data is shown when the device is started, offering the user their comprehensive measurement history at a glance.

Newly designed e-Cuff

A newly designed semi-hard shell cuff promises easy fastening around your arm as well as comfortable blood pressure measurements for patients in a wider arm range (17 to 32cm).

Unique heart shaped cuff holder

The UA-782 is equipped with a convenient cuff holder. Simply place the cuff in the cuff holder located on the backside of the unit for storage after the measurement is completed.

A&D

A&D Medical

Tokyo JAPAN

<http://www.aandd.jp>

Technology You Can Trust

*UA782CA-ADCC-01-PR1-09410
REF*CARD-UA782