


# QUICK USB PORTABLE ADAPTOR

## Model: AD-8527

### General Specifications

- Transmits weighing data directly from an A&D balance to a USB port of a PC (no driver software is required)
- Sends weighing data to software applications such as Microsoft Excel, Word and Notepad
- Small portable size
- IP65 dust and waterproof rating with casing cover



### AD-8527 SPECIFICATIONS

Applicable Devices	Balances and scales with a RS-232C connector manufactured by A&D
Dedicated power supply	Not required. Power is supplied by the weighing instrument or PC
Operating environment	5 °C to 40 °C, 85% RH or less (no condensation)
Dimensions	55mm x 103mm x 16.5mm
Weight	Approximately 60g
Standard accessories	USB extension cable (0.5m), Three RS-232C cables (1m)

\* Specifications are subject to change without notice