

Simple, Convenient Wireless Communication

# AD-8529 Series

Enables wireless communication between an A&D balance/scale and a PC/printer

Useful in fume hoods or safety enclosures with limited space and risk of contamination

Convenient where cabling is troublesome such as inside glove boxes

Allows you to prevent substances from adhering to the paper by placing a printer in a remote area



for PC



AD-8529PC-W



for Printer



AD-8529PR-W



## ■ Features

- Maximum communication range: Approximately 10m
- How to use: Plug & play  
(For the AD-8529PC-W, the driver will be automatically installed on PCs with Windows 10)
- Power supply: No battery required (excluding some models)
- Data transmission: Computer software such as WinCT

## ■ PC Application

- Connect directly to the USB port of PC
- Install the driver and WinCT software from the attached CD-Rom

## ■ Printer Application

- Multi-printer AD-8217 : Connect directly to the back of the printer
- Mini-printer AD-8216: Connection cable (AX-KO5363-30) is required

## ■ Balance / Scale Application

- D-Sub 9Pin Balance (BM, GX-A, GF-A, EK-I, etc.)  
Connect directly to the balance
- D-Sub 25Pin Balance (GR, GX, GF, GP, etc.)  
Use D-Sub 25Pin to D-Sub 9Pin cable (AX-KO5363-30)
- DIN Balance (HC-I, HW-KGV, HV-KGV, FG)  
Use DIN to D-Sub 9Pin cable(AX-KO5543-100)

※ Balance that cannot supply power to the Bluetooth converter from RS-232C

Use the AC adaptor (AX-TB253) to supply power

## ■ Options & Accessories

- AX-KO5363-30: Connection cable for a weighing device of D-sub 25pin (30cm)
- AX-KO5543-100: Connection cable for a weighing device of DIN7 Pin(1m)
- AX-TB253: AC Adaptor 100-240V



AD-8529PC-W

AD-8529PR-W

AD-8529PC/PR-W

AX-KO-5363-30

AX-KO5543-100

AX-TB253



### A&D Company, Limited

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013 JAPAN  
Telephone: [81](3) 5391-6132 Fax: [81](3) 5391-1566  
<http://www.aandd.jp>

### A&D ENGINEERING, INC.

1756 Automation Parkway, San Jose, CA 95131 U.S.A.  
Telephone: [1](408) 263-5333 Fax: [1](408) 263-0119

### A&D Australasia Pty Ltd.

32 Dew Street, Thebarton, South Australia 5031 AUSTRALIA  
Telephone: [61](8) 8301-8100 Fax: [61](8) 8352-7409

### A&D INSTRUMENTS LTD.

Unit 24/26 Blacklands Way Abingdon Business Park,  
Abingdon, Oxon OX14 1DY UNITED KINGDOM  
Telephone: [44](1235) 550420 Fax: [44](1235) 550485

### A&D KOREA Limited

817, Manhattan Bldg., 33, Gukjegeumyung-ro 6-gil,  
Yeongdeungpo-gu, Seoul, 07331, KOREA  
Telephone: [82](2) 780-4101 Fax: [82](2) 782-4280

### A&D RUS CO., LTD.

Vereyskaya str.17, Moscow, 121357 RUSSIA  
Telephone: [7] (495) 937-33-44 Fax: [7] (495) 937-55-66

### A&D Instruments India Private Limited

509 Udyog Vihar Phase V  
Gurgaon-122 016, Haryana, INDIA  
Telephone: [91](124) 471-5555 Fax: [91](124) 471-5599