

# AD-8127

## MATRIX PRINTER





### Features

- Compact and lightweight.
- Keyboard data entry.
- Highly reliable and durable.
- Compatible with all A&D scales and balances, and a wide range of other electronic balances, indicators.
- Statistical calculation and clock functions.
- Prints by keys on balances, indicators, etc. or the printer key on AD-8127, Automatic Printing Mode and Interval Printing Mode (within the range of time interval from 5 seconds to 60 minutes printing).
- Simple printer configuration.
- Provides External INPUT Signal to print by pressing an external switch.
- Self-diagnostic function.
- Customization: special communication formats available upon request.

Item	Specification
Printing method	8 Pin Serial Impact dot matrix system
Number of dot per line	240 dot
Printing speed	1.6 line/sec
Columns number	English standard 30 Columns
Font size	English: 8 x 14, Korean: 16 x 14
Support language	English, Hangul support (Optional feature by inside mode establishment)
LCD	8 x 2 Line (Back Light Function)
Printing function	Statistics arithmetic function, watch function
Cartridge	EPSON ERC-09
Printer Paper	57mm(w) x 60mm(dia) Roll-Paper
Printing width	48mm
Interface	- SERIAL RS-232C, Current Loop - RS-485(OPTION) BAUDRATE = 300/600/1200/2400/4800/9600/19200bps PARITY = None/ Even/ Odd DATA WIDTH = 8/7 bit STOP BIT = 1/2 STOP 2 External Input Port
Power source	AC 100 - 240V (AC Free Voltage)
Working conditions	-15° - 50°
Storage condition	-20° - 70°
Certificate	KC Certification, CE Certification
Size (mm)	174.3(W) x 141.3(D) x 86.9(H)

### Product size



**A&D**  
A&D Australasia

"Excellence in Measurement"

ACN 007 556 809

Head Office  
32 Dew Street  
THEBARTON South Australia 5031  
Telephone (08) 8301 8100  
Facsimile (08) 8352 7409

Victorian Office  
39 Bakehouse Road  
KENSINGTON Victoria 3031  
Telephone (03) 9372 1522  
Facsimile (03) 9372 1193

Email: [sales@andaustralasia.com.au](mailto:sales@andaustralasia.com.au)  
Internet page: [www.andaustralasia.com.au](http://www.andaustralasia.com.au)

New South Wales Office  
Unit 4, 14 Abbott Road  
SEVEN HILLS New South Wales 2147  
Telephone (02) 9674 5466  
Facsimile (02) 9674 2544

Your Dealer